


# uCentum™

comprehensive posterior system

# uBase™

iliosacral fixation system




Safety for user  
and patient

# uCentum™

## System


- Posterior stabilization for the thoracic and lumbar spine
- Complete system for open and percutaneous applications
- Rods, Ø 6.0 mm, straight and curved in titanium alloy and partly CoCr
- Polyaxial and monoaxial screws
- Pre-fixation screws and longarm screws
- Can be combined with Iliosacral screws uBase™



Cannulated and perforated screws (exception Ø 4.5 mm)



Sleeve for percutaneous application



Offset connector

Low screw profile

High degree of flexibility due to various connector lengths and additional connection options

# uBase™

## System

- Stabilization of the sacrum and ilium in combination with uCentum™
- Purely sacral, sacroiliac or purely iliac screw placement
- Extensive preparation instruments
- Cannulated screws for secure screw placement
- Navigation option

## Advantages

- High stability level with fixed-angle screws
- Low screw profile
- Safety for user and patient due to cannulated screws and navigation option



Purely sacral screw placement



Sacroiliac screw placement



Purely iliac screw placement



## Advantages

- One system for different treatment options
- Easy to use due to one basic instrument set
- Safety for user and patient due to cannulated screws
- Cement augmentation through the screws thanks to lateral perforation



uCentum™

# Components

Implants	Product number
Screw, polyaxial, Ø 4.5 mm, length 25–50 mm	CS 3802-045-025 – CS 3802-045-050
Screw, polyaxial, Ø 5.5 mm, length 25–55 mm, perforated	CS 3802-055-025 – CS 3802-055-055
Screw, polyaxial, Ø 6.5 mm, length 30–60 mm, perforated	CS 3802-065-030 – CS 3802-065-060
Screw, polyaxial, Ø 7.5 mm, length 35–60 mm, perforated	CS 3802-075-035 – CS 3802-075-060
Screw, polyaxial, Ø 8.5 mm, length 35–60 mm, perforated	CS 3802-085-035 – CS 3802-085-060
Screw, polyaxial, Ø 10 mm, length 45–60 mm, perforated	CS 3802-100-045 – CS 3802-100-060
Screw, monoaxial, Ø 4.5 mm, length 25–50 mm	CS 3803-045-025 – CS 3803-045-050
Screw, monoaxial, Ø 5.5 mm, length 25–55 mm, perforated	CS 3803-055-025 – CS 3803-055-055
Screw, monoaxial, Ø 6.5 mm, length 30–60 mm, perforated	CS 3803-065-030 – CS 3803-065-060
Screw, monoaxial, Ø 7.5 mm, length 35–60 mm, perforated	CS 3803-075-035 – CS 3803-075-060
Screw, monoaxial, Ø 8.5 mm, length 35–60 mm, perforated	CS 3803-085-035 – CS 3803-085-060
Screw, pre-fixation, Ø 4.5 mm, length 25–50 mm	CS 3804-045-025 – CS 3804-045-050
Screw, pre-fixation, Ø 5.5 mm, length 25–55 mm, perforated	CS 3804-055-025 – CS 3804-055-055
Screw, pre-fixation, Ø 6.5 mm, length 30–60 mm, perforated	CS 3804-065-030 – CS 3804-065-060
Screw, pre-fixation, Ø 7.5 mm, length 35–60 mm, perforated	CS 3804-075-035 – CS 3804-075-060
Screw, pre-fixation, Ø 8.5 mm, length 35–60 mm, perforated	CS 3804-085-035 – CS 3804-085-060
Screw, pre-fixation, Ø 10 mm, length 45–60 mm, perforated	CS 3804-100-045 – CS 3804-100-060
Longarm screw, pre-fixation, Ø 5.5 mm, length 25–55 mm, perforated	CS 3805-055-025 – CS 3805-055-055
Longarm screw, pre-fixation, Ø 6.5 mm, length 30–60 mm, perforated	CS 3805-065-030 – CS 3805-065-060
Longarm screw, pre-fixation, Ø 7.5 mm, length 35–60 mm, perforated	CS 3805-075-035 – CS 3805-075-060
Rod, Ø 6.0 mm, straight, length 30–500 mm, standard approach	CS 3807-030 – CS 3807-500
Rod, Ø 6.0 mm, straight, length 30–500 mm, percutaneous approach	CS 3808-030 – CS 3808-500
Rod, Ø 6.0 mm, curved, length 30–80 mm, percutaneous approach	CS 3809-030 – CS 3809-080
Rod, CoCr, Ø 6.0 mm, straight, length 400 mm, hex end, standard approach	CS 3812-400
<b>Combination with uBase™</b>	<b>Product number</b>
Iliosacral screw, Ø 7.5 mm, length 40–80 mm	CS 7170-075-040 – CS 7170-075-080
Iliosacral screw, Ø 8.5 mm, length 60–120 mm	CS 7170-085-060 – CS 7170-085-120
Iliosacral screw, Ø 10.0 mm, length 80–140 mm	CS 7170-100-080 – CS 7170-100-140

Over a Century  
of Innovation



ulrich GmbH & Co. KG | Buchbrunnweg 12 | 89081 Ulm | Germany  
Phone: +49 (0)731 9654-0 | Fax: +49 (0)731 9654-2705  
spine@ulrichmedical.com | www.ulrichmedical.com